

Core Motion Capture Workflow (OpenCap – Treadmill Setup)

What You Need


- **1 Laptop** with stable WiFi
 - **2 iOS devices** (2018 or newer, on the same WiFi)
 - **OpenCap app** installed on both mobile devices
 - **Checkerboard** (flat or perpendicular to the floor, to avoid calibration errors)
 - **Credentials file for authentication in opencap.ai** (mocap_cred.txt)
-

Setup & Camera Positioning

1. Place cameras **on both sides of the treadmill**, each at **~45° angles** from the subject's midline (see sketch in appendix).
 2. Position cameras at **hip to shoulder height**, ensuring the **entire body is in view**.
 3. Ensure subjects wear **fitted clothing** (test edge cases with coats, skirts, hijabs if possible).
-

Authentication & Session Setup

1. Log in: app.opencap.ai/login using the credentials in file.



2. Create a **NEW SESSION** and follow prompts.
 3. Mount and align both mobile devices using the QR codes and app responses.
-

Calibration & Consent

4. Run **data assessment + calibration** for each subject.
 - a. Record approximate **height and weight**.
 - b. Obtain **consent** (data is processed in the cloud, videos may be de-identified).
-

Recording Trials

5. Record **short treadmill trials** per subject (trial length: 30 to 60 seconds):
 - a. At least **2 normal walking** trials.

- b. **2–5 varied trials** (different speeds or intentional gait anomalies like stiff knee, hip misalignment, etc.).
 6. Keep **subject fully visible** at all times.
 - a. Avoid multiple people entering the frame (causes inverse kinematics errors down the line processing).
-

✓ Key Success Factors for Recording

- Reliable WiFi + properly mounted iOS devices.
 - Correct camera placement (angles, height, and full-body visibility).
 - Calibration with checkerboard is essential to avoid skeleton drift.
 - Subject consent for cloud-based data processing.
 - Multiple trial recordings for robust dataset, including normal and edge-case movements.
-

Setup Sketch:

